
RIKTA®

RIKTA®-VET

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EN**ATTENTION!**

- When purchasing RIKTA®-VET device check the following:
 - the complete set of delivery (p. 26);
 - compliance of serial numbers in this Operating Manual (p. 75) with the numbers on the device;
 - availability of sale date and stamp of the seller in the Warranty Card (p. 75).
- **Do not switch on the device if it was kept at temperatures below 0 °C!** The device should be kept in conditions of room temperature for at least 2 hours.
- Before operating this device, read this Manual carefully and keep it during the entire useful life of the device.

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Dear Customer!

Congratulations on your purchase of the RIKTA®-VET device!

We express our confidence that the device will help your pets to keep their bones green for many years to come!

INVISIBLE LASER RADIATION!
DO NOT STARE INTO THE BEAM AND DO NOT MAKE
DIRECT OBSERVATIONS USING OPTICAL INSTRUMENTS!
CLASS 1M LASER EQUIPMENT

TABLE OF CONTENTS

RU 3

EN

- GENERAL INFORMATION 22
- SAFETY MEASURES 24
- DELIVERY SET 26
- TECHNICAL PARAMETERS OF RIKTA®-VET DEVICE 27
- DEVICE 28
- PREPARATION FOR OPERATION 30
 - Device turning on 30
 - Device serviceability check (testing) 31
- HOLDING TREATMENT SESSIONS 32
- DEVICE TURNING OFF 33
- BATTERY CONDITION MONITORING 34
 - Battery recharging 34
- MAINTENANCE 35
- COMMON TROUBLES AND REMEDIES 36
- PACKAGING, STORAGE AND TRANSPORTATION 37
- MANUFACTURER'S GUARANTEE 38

FR 39

DE 57

- ACCEPTANCE CERTIFICATE 75
- GUARANTEE CARD 75

The multi-year research of large group of scientists and veterinarians enabled creation of a portable magnetic-infrared laser therapeutic device for veterinary use RIKTA[®]-VET (hereafter, the Device).

The device is produced by Russian enterprise “MILTA – PKP GIT” JSC which quality control system complies with requirements of international standards.

RIKTA[®]-VET Device meets requirements of European and international electrical safety and electromagnetic compatibility standards.

The Device is certified in the Russian Federation.

Certificate of Conformity No. PCCC RU.HA34.H01233.

GENERAL INFORMATION

- The Operating Manual (hereafter, the **Manual**) is a document that certifies the manufacturer-guaranteed basic parameters and technical characteristics of the magnetic-infrared laser therapeutic device RIKTA[®]-VET. The Manual is intended for familiarization with the operation rules and is a guide for maintenance, transportation, and storage.

**The Manual shall be kept safe during
the entire useful life of the device.**

- The treatment procedure methods are given in the “Guidance manual on the use of the quantum therapy device RIKTA[®]-VET” (hereafter, the **Guidance Manual**) included in delivery of each Device.
- The manufacturer reserves the right to make changes and amendments to this Manual and the technical documentation at any time without notice.
- The manufacturer reserves the right to make changes in the complete set and the design of the Device that do not affect its functioning and basic technical parameters.

PURPOSE

The device is designed for non-drug, pain-free and non-invasive (without skin penetration) treatment of a wide range of animal diseases. Combination of the laser therapy using RIKTA®-VET Device with the pharmaceutical treatment methods improves the efficiency of drug administration.

The Device provides simultaneous deep penetration of the near infrared laser radiation into the tissue in combination with the constant magnetic field, as well as light fluxes of infrared, red and blue band.

The treatment can be given both by a contact method (the emitting window of the Device comes in contact with skin) and by non-contact method (exposure at a distance of 0.5–2 cm from skin). In case of big sizes of the affected area, the treatment is given by scanning, i.e. uniform motion of the Device at a speed of 0.5–1.5 cm/s.

The Device is designed for use in veterinary clinics, as well as for treatment of animals in the field and domiciliary.

SAFETY MEASURES

The Device may be operated only after familiarization with this Manual and the Guidance Manual.

DO NOT BREAK THE SEAL

- Never open the Device and never try to disassemble or modify its internal components. The Device has no components serviceable by the user. In case of malfunction, immediately stop using the Device and resort to the manufacturer or the service center.

PLUG AND POWER CORD

- Do not plug the battery charger into mains with wet hands to avoid electric shock.
- Do not pull the battery charger cord when disconnecting the power cord, make sure to hold the plug when disconnecting the power cord.

WHEN MALFUNCTIONS ARE OBSERVED

- In case of wear or damage of the cord or battery charger plug, as well as in case of strange odor and smoke, immediately turn off the Device, disconnect the power cord from mains and resort to the manufacturer or the service center.

EXTERNAL EFFECTS

- It is prohibited to spray flammable substances in the vicinity of the Device.
- Avoid exposure of the Device to moisture and excessive heat.
- The device complies with ISO 60601-1-2 standard for electromagnetic compatibility. However, as a precaution it is not recommended to use the Device nearby the devices sensitive to electromagnetic radiation (e.g. mobile phones, medical equipment, computers, etc.).

WHEN USING THE DEVICE

- The personnel should wear goggles (optionally). Single treatments can be given without goggles provided that the emitting window of the Device tightly adjoins the animal body.
- The Device falls within Class 1M laser equipment.


When using the Device, it is recommended to avoid eyes exposure to direct or reflected radiation. Do not stare into the beam and do not make direct observations using optical instruments.

- Start a treatment procedure only when the emitting window is targeted at the affected area.

DEVICE CARE

- Before and after a treatment procedure, make sure to wipe the Device with a cloth of gauze or coarse calico soaked in a treatment tool disinfectant based on 3% hydrogen peroxide solution. Avoid Device contact with alcohol-containing media!
- The Device disinfection may be carried out only when the Device is turned off and disconnected from mains.

OPERATION CONDITIONS

The Device is designed for operation at a temperature from +10 to +35 °C, atmospheric pressure from 630 to 800 mm Hg, and relative humidity not more than 80%.

EN

The Device and its components should be disposed in accordance with general applicable standards and local environment protection rules.

DELIVERY SET

The complete set of Device delivery is specified in Table 1.

Table 1

No.	Description	Qty
1	RIKTA®-VET device	1 pc.
2	Operating manual	1 pc.
3	Guidance Manual	1 pc.
4	Battery charger c/w power cord	1 pc.
5	Packing bag	1 pc.

Optionally: a set of light-guide nozzles KON-1 adding functionality of the Device.

TECHNICAL PARAMETERS*

OF RIKTA® –VET DEVICE

Emission wavelength, μm :

- laser 0,800...0,910
- broadband infrared..... 0.860...0.960
- visible red radiation 0.600...0.700
- visible blue radiation..... 0.465...0.475

Pulse recurrence frequency, Hz:

- red band2
- blue band.....2
- infrared and laser radiation:
 - constant.....4; 64; 512; 4096
 - variable (VAR)..... 250÷1

Treatment procedure duration, min2; 5; 10

Pulse power of laser radiation, W 21 ± 1

Average power of broadband infrared radiation, mW..... 60 ± 30

Magnetic field induction, mT 35 ± 10

Emitting window area, cm^2 $4 \pm 0,4$

Power supply from internal Li-Po battery:

- voltage, V 3.7
- ampere-hour capacity, mA•h1800

Continuous operation time with fully charged battery, hour, min8

Battery full charging time, hour, max.....3

Output voltage of AC adapter for battery charging, V 13.5

AC adapter input voltage 110...240 V / 50 Hz

Laser safety according to IEC 60825-1 Class 1M

Electrical safety according to IEC 60601-1 BF type


Overall dimensions, mm 203 x 64 x 70

Net weight, g, max (w/o adapter)..... 250

Service life of the Device, year5

* For reference

DEVICE

RIKTA®-VET Device

The Device body is made of impact-resistant and environmentally friendly plastic.

The Device runs on the built-in battery.

- 1 – laser radiation frequency select button (**FREQUENCY**) or red light turn on / off button (**RED ON/OFF**);
- 2 – laser radiation frequency indicators: 4 Hz, 64 Hz, 512 Hz, 4096 Hz, VAR – variable frequency (250...1 Hz);
- 3 – battery charge indicator (**Battery**);
- 4 – treatment session time indicators: 2 min., 5 min., 10 min.;
- 5 – treatment session time setting button (**TIME**) or blue light turn on / off button (**BLUE ON/OFF**);
- 6 – laser radiation presence sign;
- 7 – Device turn on / off button or treatment session start / stop button (**START/STOP**);
- 8 – thread for light-guide nozzles;
- 9 – emitting window;
- 10 – battery charger in jack;
- 11 – Device nameplate:


- T1 – manufacturer's logo;
- T2 – certification mark of the RF Federal Agency on Technical Regulating and Metrology;
- T3 – laser radiation presence sign;
- T4 – laser safety class;
- T5 – electrical safety mark, BF type;
- T6 – No. of technical specifications which the Device conforms to;
- T7 – output parameters of AC adapter for battery charging;
- T8 – positive polarity mark;
- T9 – serial No. of the Device;
- T10 – manufacturer.


PREPARATION FOR OPERATION

After transportation at temperatures below zero, the Device should be kept at a room temperature of $+20 \pm 5 \text{ }^{\circ}\text{C}$ during 2 hours.

DEVICE TURNING ON


Turn on the Device by pressing and holding in position **START/STOP** button until an audible signal is heard.

Holding
treatment
sessions
page 32

If after turning on the Device **only “64 Hz” and “5 min” indicators light up**, it is possible to commence the treatment procedure (see “Holding Treatment Session”, page 32).

Service-
ability
check
page 31

If the Device had been off for more than 45 minutes then when the Device is turned on **all frequency indicators [2] and all time indicators [4]** will light up. At that **FREQUENCY** and **TIME** buttons would not respond to pressing: the Device is standing by waiting on the serviceability check (see “Device Serviceability Check”, page 31).

If after the Device is turned on or if during the Device operation **“Battery” [3] indicator starts blinking red**, it means the low battery, so for further normal operation of the Device it is necessary to recharge the battery (see “Battery Recharging”, page 34).


DEVICE SERVICEABILITY CHECK (TESTING)

If the Device had been off for more than 45 minutes then when the Device is turned on, the serviceability check standby mode will be activated automatically. At that all frequency [2] indicators and all time [4] indicators will light up. **FREQUENCY** and **TIME** buttons would not respond to pressing.

Make sure that all frequency [2] and time [4] indicators are on but do not blink!

If the indicators blink the serviceability check can not be carried out. As a rule, it happens if the emitting window [9] of the Device is targeted at a light source. Position the Device so that the emitting window [9] is not targeted at a light source. When the indicators stop blinking it is possible to proceed to the Device serviceability check.

For the Device serviceability check:

If after the audible signal only “64 Hz” and “5 min” indicators are on **the Device is serviceable and ready for operation** (see “Holding Treatment Session”, page 32).

If all frequency [2] and time [4] indicators blink and intermittent audible signal is heard it means that **the Device is faulty**. Immediately turn off the Device and resort to the manufacturer or the service center.

HOLDING TREATMENT SESSIONS


The sequence of the session, the exposed areas, the laser radiation frequency and the exposure time are specified in the Guidance Manual or in an individual order of a veterinarian.

1

In accordance with recommendations set forth in the Guidance Manual or in the order of a veterinarian:

- sequentially pressing **FREQUENCY** button, set the required laser radiation frequency: **4 Hz, 64 Hz, 512 Hz, 4096 Hz, VAR** (variable frequency);

- sequentially pressing **TIME** button, set the exposure time: **2 minutes, 5 minutes, 10 minutes.**

2

Target the emitting window of the Device [9] at the specified exposure area.

Press **START/STOP** button.

An audible signal will be heard:

treatment of the area has been started.

From the exposure start and until the exposure completion, **FREQUENCY** and **TIME** buttons will operate, accordingly, in the mode:

“**RED ON/OFF**” (turn on/off red light radiation)

“**BLUE ON/OFF**” (turn on/off blue light radiation).

**At the moment of the treatment procedure launching,
red light radiation is on,
blue light radiation is off.**

During a treatment session, the indicator of selected frequency blinks, and the time indicator is permanently on.

Upon the expiry of the set time, a triple audible signal will be heard: **treatment of the area has been completed.**

**A treatment procedure can be interrupted at any time
by brief pressing on START/STOP button.**

DEVICE TURNING OFF

Upon completion of the Device use make sure to turn it off: press and hold **START/STOP** button for 3–5 seconds

All indicators will go out, and the Device will be turned off.

**If in the “on” condition the Device is not used
(no buttons are pressed, no treatment is carried out) for 40 seconds,
the Device will be turned off automatically.**

BATTERY CONDITION MONITORING

After turning on the Device and during Device use, the battery condition is monitored automatically.

“**Battery**” indicator [3] shows the state of battery charge:


- If the battery voltage is normal the “**Battery**” indicator [3] is off.
- If “**Battery**” indicator [3] blinks red the stored energy of the battery will be consumed within max. 10 minutes of Device operation, so the battery should be recharged!

In case of low battery the Device will be turned off.


BATTERY RECHARGING

In order to charge the battery:

1. Plug the battery charger connecting cable in the jack of the Device

2. Plug the power cord in the battery charger jack

3. Insert the power cable plug in the main socket

During the charging process “**Battery**” indicator [3] lights red.

Once the battery is charged (“**Battery**” [3] indicator will light green), disconnect the power cable from (mains and remove the connecting cable from the Device jack.

MAINTENANCE

- In order to ensure reliable operation of the Device it is necessary to perform the proper maintenance as specified in Table 2 hereof.
- During the maintenance, make sure to take safety measures specified in section “SAFETY MEASURES”, page 24.
- If the Device is found non-conforming to the technical requirements specified in Table 2, it may not be operated any more and is subject to repair.

Table 2

Maintenance Frequency	Work Scope	Technical Requirements
Before each session	Visual inspection not using any special tools	No mechanical damages on the housing of Device and battery charger
	Check the following:	
	1. State of battery charge by Battery indicator [3]	Battery indicator [3] should be off
	2. Availability of radiation at the emitting window output [9]	Once the emitting window [9] of the Device is targeted at the exposure area and START/STOP button is pressed, visible red light radiation should be observed
3. Disinfection (if required)	According to recommendations given in section “SAFETY MEASURES”, page 24	

COMMON TROUBLES AND REMEDIES

The list of common troubles, probable causes and the remedies are specified in Table 3.

Table 3

Trouble	Probable cause	Remedy
<p>After START/STOP button is pressed to turn on the Device, indicators [2] and [4] do not light up.</p>	<p>Low battery or faulty Device</p>	<p>Recharge the battery (see section “BATTERY CONDITION MONITORING”, page 34).</p> <p>If after battery charging and turning on the Device, [2], [4] indicators do not light up, the repair works shall be carried out by the manufacturer.</p>
<p>When checking the Device serviceability (see section “Device Serviceability Check”, page 31) all [2], [4] indicators blink, an intermittent audible signal is heard</p>	<p>One or more types of radiation are missing</p>	<p>The repair works shall be carried out by the manufacturer</p>

PACKAGING, STORAGE AND TRANSPORTATION

- Each Device is delivered together with all accessories in the package of the manufacturer.
- The Device shall be stored in the package of the manufacturer in a heated space at a temperature of +10...+35 °C, air relative humidity max. 80 % (at a temperature of +25 °C). The storage space shall be free from corrosive dust, mercury vapors, acids, and alkalis.
- The Device set may be transported by any enclosed transport in the package of the manufacturer. It is necessary to ensure a stable position of containers during transportation.

MANUFACTURER'S GUARANTEE

- The manufacturer guarantees the conformity of the Device to requirements of technical specification TY 9440-004-18537172-2002 during 24 months from the sale date of the Device to the consumer (subject to observation of operation, transportation and storage conditions by the consumer). During the specified period, the manufacturer carries out free-of-charge repair of the Device.
- The manufacturer bears no responsibility and carries out no free repair when troubles occurred through the fault of the consumer (for instance, in result of careless handling, transportation, failure to observe the instructions for use, as well as in case of the broken seal).
- The Device guarantee does not apply in case of:
 - unavailability of the Manual;
 - mechanical damages, including damages occurred during transportation;
 - failure of the Device because of foreign bodies or liquids found inside the Device;
 - broken seal.
- The guarantee shall be void if during the warranty period the repair of faulty Device was carried out by unauthorized persons and/or organizations.
- Upon expiry of the guarantee period, the repair of the Device should be carried out at the expense of the consumer.
- If the guarantee card bears no sale date and seller's stamp, the guarantee period starts from the Device manufacture date.
- As far as repair and operation issues are concerned, the consumer should contact the manufacturer and service centers specified on the web-site www.rikta.ru or in separate guarantee card (if any).

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